

MANO XXL

The super large opto-electronic magnifier
with front camera



MANO XXL offers a super large display of 10" in attractive design and a modern look like a tablet with a continuous magnification of up to 32x.

Due to the auto focus and the high resolution FullHD camera, all text and pictures are being displayed in razor sharp quality and contrast enhanced. Take a snapshot and magnify it in the comfort of your individual needs, and/or change contrast. An additional front camera brings you into the position to look at yourself while hairdressing or doing

make up. Performing handwritten tasks, using distance viewing and/or mirror function are just a few more features that can be comfortably done by using the comfortable reading/writing stand.

Intuitive navigation via joystick is comfortable and makes operation easy.

Connect MANO XXL with an external monitor or TV via the integrated HDMI port and look at your photos or documents at an even higher magnification level.

Your benefits at a glance

- Super large 10" display (approx. 25,4 cm diagonal)
- continuous magnification from approx. 2x to 32x
- razor-sharp image in all magnification levels thanks to FullHD camera and autofocus
- intuitive operation via joystick, snapshot function, line- and column marker
- 2 cameras (reading- and front camera)
- integrated stand, foldable reading and writing stand
- HDMI port, USB port



read and write
comfortably



MANO XXL

order code:
MANXXLP3



10" display (approx. 25,4 cm diagonal)
snapshot function, memory of approx.
2GB, autofocus



approx. 2-32x
system resolution 1280 x 800 pixel



Lithium-Polymer-battery
operating time approx. 4.5-5.5h,
charging time approx. 5.5 h




true colors,
up to 16 selectable semi colors



weight approx. 700 g
W x H x D approx. 24,5 x 20 x 2,2 cm



scope of supply: reading and writing
stand, carrying bag, power supply with
country specific cord, HDMI cable, user
manual



folding reading
stand

second
camera



Reinecker Vision GmbH
Head office Germany
Sandwiesenstraße 19
D-64665 Alsbach-Hähnlein

Your authorized distributor:

Tel. +49 6257 9311-0

Fax +49 6257 9311-8008
international@reineckervision.de