

VEO EDU

With distance camera and split-screen function



Your benefits at a glance

- ideal for use at school and at work
- LVHD technology for high-contrast, flicker-free picture in FullHD quality
- high resolution 24" class LED flat screen (approx. 60 cm)
- movable monitor with spring compensation for comfortable reading posture and easy positioning
- intuitive 3-button operation of the document camera
- electrically adjustable distance camera with freely positionable control element
- 2 HDMI inputs for displaying external image sources (PC/smartboard/additional camera)
- USB 3 output enables the VEO image to be displayed and processed (saved/read out) on a PC or laptop
- display of up to four different image sources (2 image sources simultaneously)

The Focus camera

- FullHD camera with digital zoom
- magnification approx. 3.7x to 42x (factory setting, extended magnification range (starting at approx. 2.2x) can be set in the menu)
- razor-sharp image even at a low zoom level

The Vario distance camera

- FullHD+ distance camera with telescope stand, separately controllable, memory function for four camera positions
- high ease of use: dynamic camera movement to the desired position

VEO EDU

order code:
VEVED244



LED flat screen monitor
24" class | approx. 60 cm
(1920x1080 px)
autofocus, can be switched off



FullHD document camera
with digital zoom
electronically controllable FullHD+ distance camera with 18x optical zoom



magnification approx. 3.7 to 42x*
*Factory setting, extended magnification range can be set in the menu



true colors,
semi colors: 
user-defined semi color modes



Weight approx. 13.2 kg
Dimensions (WxHxD)
approx. 54 x 56 x 63 cm



Scope of supply: distance camera control panel with joystick, power supply unit with cable



Accessories (optional):
height adjustable table with control panel holder and cable storage

control panel for the distance camera, split screen, etc.



comfortable reading with XY table

simple 3-button operation of the document camera

motorized distance camera

height-adjustable from approx. 57 to 79 cm

USB 3 output

Connect up to 2 external image sources via HDMI

Reinecker Vision GmbH

Head office Germany

Sandwiesenstraße 19

D-64665 Alsbach-Hähnlein

Tel. +49 6257 9311-0

Fax +49 6257 9311-8008

international@reineckervision.de

Your authorized distributor: