

mezzo

The compact desktop video magnifier



Whether you choose to use the mezzo as a compact desktop video magnifier, or as a classic reading system with an XY table, the mezzo performs with crystal clear, sharp, high-resolution pictures even at the highest magnification levels.

The FullHD camera enlarges your reading material with the best picture quality, minimal picture noise and reduced reflection of glossy surfaces. You can also benefit from reduced shadows while hand-writing.

The high resolution, LED, flat screen monitor, is available in different sizes. With its ergonomic, movability, and spring form positioning, the mezzo is truly the most comfortable system available.

With its compact design and easy folding mechanism you can be assured that transporting and set up will be a breeze.

mezzo Focus 24" comes always with an integrated XY table with locking mechanism, whereas mezzo Focus 20" can be ordered with or without XY table.

Your benefits at a glance

- **LVHD**- Technology: high resolution with virtually no picture noise
- high resolution LED flat screen monitors available in two sizes: 20" and 24" class
- reduced shadows while hand-writing
- high resolution image in HD quality
- movable monitor with spring form positioning for an ergonomic reading posture
- effortless set up and folding for transportation or storage
- easy to operate three button control

The Focus camera

- FullHD camera with digital zoom
- excellent image quality even when using a low zoom level
- magnification approx. 2.1x to 45x (dependent on monitor size, options and settings)

mezzo

order codes: MEZFO200/
MEZF20X1/MEZF24X1



LED flat screen

20" class: 1600 x 900 pixel
24" class: 1920 x 1080 pixel

comfortable autofocus with one-touch
deactivation



Focus camera:
FullHD camera with digital zoom



20" class: approx. 2.6x to 30x*
24" class: approx. 3.8x to 45x*

*refers to factory settings, further magnification lev-
els to be adjusted within the internal menu



true colors,
semi colors:

individually selectable semi color
modes



20" class: weight approx. 6.2 kg
dimensions (W x H x D) approx.
unfolded 46 x 50 x 52 cm
folded 46 x 18 x 52 cm

24" class: weight approx. 10.5 kg
dimensions (W x H x D) approx.
unfolded 53 x 53 x 52 cm
folded 53 x 22 x 52 cm



scope of supply:
AC adapter, user manual



available options:
XY table, carrying bag



comfortable
reading with cross table
simple 3-button operation



space-saving
foldable

Reinecker Vision GmbH

Head office Germany

Sandwiesenstraße 19

D-64665 Alsbach-Hähnlein

Tel. +49 6257 9311-0

Fax +49 6257 9311-8008

international@reineckervision.de

Your authorized distributor: