

topolino flex 2K

The camera system with 2K resolution for more working space



The camera system topolino flex 2K offers you the greatest possible working space on your desk. The high-resolution camera (2.5x HD resolution) is attached to a three-joint swivel arm on the desk. This allows you to move the camera at different angles across the entire workspace and display the image on a monitor. All individual adjustments in terms of colors, magnification, etc. are being saved automatically when switched from reading mode to distance viewing mode.

Via an HDMI interface, topolino flex can be connected directly to a monitor or connected to a PC or laptop via USB3.0. Your external 2K monitor can be connected via HDMI interface with topolino flex 2K. Your PC or laptop can be connected via USB 3.0. Functions for your comfort, like line ruler or line masking can be activated via topolino remote control and the Viewer Software SuperNova for Window.

Your benefits at a glance

- space saving, flexible 2K reading-/ distance viewing camera with a highest resolution of up to 2048 x 1280p*
- approx. 19x lossless magnification in FullHD, magnification in reading mode of approx. 0.8x to 80x*
- compatible with Windows 10, 11, macOS und iPad
- ergonomic, individual adjustable reading position
- razor sharp display of pages larger than DIN A3
- intuitive, flexible operation: via camera, remote control or PC keyboard

*depending on monitor size and connection (USB 3.0 or HDMI)

topolino flex 2K

order code:
TOFL2K31



60 Hz camera, progressive scan, system resolution of up to 2048 x 1280 pixels in distance viewing mode
autofocus (can be switched ON/OFF)



approx. 0.8x to 80x in reading mode
*depending on monitor size and connection (USB 3.0 or HDMI)



true colors
semi colors: **F F F F F**
individually selectable semi color modes



weight approx. 4,5 kg
dimensions: (W x H)
approx. 26 x 58,5 cm,
radius of action approx.
45 to 85 cm



compatible with Windows 10 and 11,
macOS, iPad (with USB-C port)
connections: HDMI port, USB 3.0 port



scope of supply: remote control, power supply, USB 3.0 cable



available options: mobile XY table, XY table with friction brake;
PC, notebook, monitor,
PC magnification software

easy
operation



Reinecker Vision GmbH
Head office Germany
Sandwiesenstraße 19
D-64665 Alsbach-Hähnlein

Your authorized distributor:

Tel. +49 6257 9311-0

Fax +49 6257 9311-8008
international@reineckervision.de