

topolino pad

Turns your iPad into a screen reader with a distance camera and read-aloud function



The topolino pad is the perfect companion for your iPad and offers excellent magnification for people with visual impairments. Whether in class, in the office or at home – with the topolino pad, you have everything in crystal-clear focus.

Thanks to the intuitive controls, you can effortlessly switch between near and far vision. Ideal for lectures, presentations, note-taking and much more. Use the

convenient split screen function to display both camera images at the same time.


The integrated read-aloud function supports over 20 languages.

With a battery life of over 5 hours, you are also flexible when you are on the go. The foldable stand allows for safe transport together with your iPad and can be securely transported using the carrying handle.

Your benefits at a glance

- thanks to a 7x optical zoom with up to 32x magnification, even the smallest details can be displayed sharply
- the full HD distance camera with 60 fps ensures a smooth and steady image, even with fast movements
- the compact device is easy to transport and can be set up in no time
- text-to-speech function available in more than 20 languages

topolino pad		order code: TOPOPAD1
	7x optical zoom, magnification approx. 1x to 32x	
	Full HD distance camera with 60 fps autofocus	
	true colors, semi colors 4 viewing options: split screen; distance view; close-up view; selfie	
	Read-aloud functions: over 20 languages included	
	Lithium-ion battery Operating time of up to approx. 5h Charging time approx. 4.5 h	
	System requirements: iPads with iPadOS16 or higher	
	Weight approx. 2.4 kg (without tablet)	
	Dimensions (W x H x D, folded) approx. 31.2 x 5.5 x 30.9 cm Height (set up): approx. 27.8 cm	
	scope of supply: carry case, power adapter with charging cable, con- nection cable for tablets, microfiber cloth, magnet holders for tablets	



read aloud
function



Reinecker Vision GmbH
Head office Germany
Sandwiesenstraße 19
D-64665 Alsbach-Hähnlein

Your authorized distributor:

Tel. +49 6257 9311-0

Fax +49 6257 9311-8008
international@reineckervision.de