

topolino tablet

The multimedia camera system that magnifies and reads aloud



topolino tablet is a compact multimedia camera system that supports pupils, students and professionals to easily fulfill their tasks. It's a 13" Windows tablet with a pre-installed magnifying- and text recognition software, that comes with a foldable stand on which a swivel-tilt distance camera can be mounted. Writing, reading text, viewing pictures, looking at objects, distance viewing, etc. can be done effortlessly and with fun. The split-screen function offers display of two different images simultaneously on the screen, e.g. distance camera and

internal camera. Magnification level and details to be displayed can be adjusted by swiping on the touch screen. Due to the integrated text recognition with read-aloud-function you can take a picture of a text and double click on the desired part of the text that you want to be read aloud.

Based on a Windows operating system, topolino tablet can be upgraded to a fully-fledged mobile workstation with a mouse and keyboard, as well as all common office applications and assistance software such as SuperNova.

Your benefits at a glance

- Windows tablet with a large 13" (approx. 33 cm) touch display on an inclinable stand
- magnification in real time and text recognition of a full DIN A4 page, split screen function
- reading camera with integrated LED illumination; distance viewing camera
- easy to fold and carry with a practical handle
- compatible with common assistive programs (e.g. SuperNova) and Office applications for Windows
- optionally available: second display for true split-screen function;
- can also be used as a notebook with large-letter keyboard and mouse

topolino tablet

order code:
TOPOTAB4



display: 13" (approx. 33 cm)
touch display on an inclinable stand



zoom: min. 1x – 60x
in document viewing mode



close-up reading camera with LED
illumination, distance viewing camera,
autofocus



Li-Ion battery, operating time up to
approx. 16 h in laptop mode



true colors, high contrast, grey scale,
semi colors:



Windows tablet, Windows 11 pro,
topolino tablet software
USB port, headphone connection



voice output: German, English
and many more languages



weight: approx. 3.5 kg

Dimensions: (w x h x d):
approx. 42 x 28 x 28 cm
folded: approx. 42 x 4 x 29 cm



scope of supply: USB-Distance
viewing camera, carrying bag,
charging cable, software full
version topolino tablet, software
trial version Dolphin SuperNova



available options: second display,
Bluetooth mouse and large characters
key board, PC magnification software,
e.g. SuperNova, Microsoft Office



Reinecker Vision GmbH

Head office Germany

Sandwiesenstraße 19

D-64665 Alsbach-Hähnlein

Tel. +49 6257 9311-0

Fax +49 6257 9311-8008

international@reineckervision.de

Your authorized distributor: