

MANO PRO

The innovative video magnifier with smartphone connection



Whether at school, at work, or at home, MANO PRO offers maximum flexibility in every situation with its split-screen function, freely movable distance camera, and the option to display an external device.

MANO PRO has a document camera and a powerful distance camera. Unlike other video magnifiers, MANO PRO can

also magnify content from your smartphone or tablet and display it in various color combinations, just like a conventional video magnifier. Android phones and tablets can even be operated directly on the screen using an app.

MANO PRO not only supports reading, but also offers text recognition and read-aloud functions in various languages.

Your benefits at a glance

- 3-in-1 functionality: document camera, distance camera, and device mirror function
- 17.3" touchscreen with matte finish
- easy connection to mobile phones or tablets for magnified display of digital content
- control Android devices directly via the MANO PRO touchscreen
- rotatable screen for horizontal and vertical viewing
- full-page mode: enlarge and navigate through documents without having to move the page
- text recognition and read-aloud function
- compact, lightweight and portable design with long battery life

MANO PRO order code: MANOPRO1	
	17,3" touch display with approx. 43.9 cm diagonal
	Magnification: document camera: approx. 1.8 to 85x distance camera up to approx. 30x
	system resolution 1920 x 1080 pixel 3 cameras for close-up/long-distance viewing in UltraHD 4K resolution (3840 x 2160 pixel)
	true colors, up to 16 selectable semi colors
	OCR and read-aloud functions: over 15 languages with different voices included
	USB-C port for smartphone/tablet, HDMI input for PC/laptop, SD card slot/USB-A port for data transfer
	replaceable lithium-ion battery operating time up to approx. 5 h (depending on type of application) charging time approx. 3 h
	weight approx. 3,6 kg
	dimensions (WxHxD) approx. 46,8x40x11,1 cm
	Scope of supply: USB-C cable, Lightning to USB-C cable, power adapter, microfiber cloth

adjustable 4K
distance camera



Smartphone display
and control
(Android only)



Reinecker Vision GmbH
Head office Germany
Sandwiesenstraße 19
D-64665 Alsbach-Hähnlein

Tel. +49 6257 9311-0

Fax +49 6257 9311-8008
international@reineckervision.de

Your authorized distributor: