

topolino smart 2K

The camera system with 2K resolution for flexible use



topolino smart 2K supports pupils, trainees, students and professionals in all tasks based on reading and thus secures integration. 2K resolution in distance viewing mode in combination with a high grade of mobility – is the ambition of this product development: moving to different classrooms, offices, etc. – with topolino smart 2K you have the right product to support you.

topolino smart 2K is being installed within seconds. All individual adjustments made are being saved when switching from distance viewing mode to reading mode and will be recalled. The dimming light is helpful when reading at different locations with different light conditions to eliminate fatigue of your eyes.

In combination with a notebook or tablet, images can be saved and videos recorded.

Your benefits at a glance

- space saving and light weight reading- / distance viewing camera system at highest resolution of up to 2048 x 1280p* (in distance viewing mode); approx. 22x lossless magnification in fullHD, magnification in reading mode approx. 1.25x to 100x*
- compatible with Windows 7, 8, 10, macOS
- dimmable illumination, battery operation,
- reproduction of a complete DIN A4 page
- automatic saving of individual adjustments in reading- / and distance viewing mode
- Windows software split screen for simultaneous display of PC image and images of reading- / distance viewing mode on a monitor
- compatible with SuperNova and Zoomtext, HDMI- and USB 3.0 interface

topolino smart 2K



camera: 2K 60 Hz camera,
progressive scan,
system resolution up to 2048 x 1280
pixels* in distance viewing mode
comfort auto-focus to be switched
ON and OFF



zoom: approx. 1.25 x to 100 x



Li-Ion battery: operating time up to
approx. 3 h, w
charging time approx. 2.5 h



true colors and artificial colors to be
adjusted via software:
F F F F F



weight: approx. 2.3 kg
Dimensions: (w x h x d) approx.:
19 x 45 x 45 cm
folded approx.: 19 x 10 x 36 cm



operating system: Windows, macOS
connections: HDMI-port, USB 3.0 port



scope of supply: carrying bag, power
supply, USB 3.0 cable, user manual



available options: mobile XY table,
comfort XY table with 2 friction brakes,
notebook, monitor, HDMI-switch (2 x
input/1 x output) with foot pedal,
PC magnification software



Reinecker Vision GmbH
Head office Germany
Sandwiesenstraße 19
D-64665 Alsbach-Hähnlein

Tel. +49 6257 9311-0

Fax +49 6257 9311-8008
international@reineckervision.de

Your authorized distributor: