

VOXbox

The modular reading system with voice



We're here for you

We at Reinecker Vision are happy and proud to keep our traditional, independent, family owned standing since 1971 and always put our eye on intelligent and innovative product developments, all made in Germany.

With all of our product innovations that we export to more than 50 countries, we focus first on the demands of our low vision individuals.

Whether in the home, educational or professional environment we design products to help our customers reach their full potential.

Your benefits:

- more than 45 years of experience in low vision solutions
- short decision making progress
- individual support
- personalized customer service

Reinecker Vision

Everything from one source

If your vision declines considerably at an older age or as a result of an eye disease, a wide range of assistive devices can help you manage your everyday life. We are happy to help you put together individual solutions to improve your daily life.

Our product range for your well-being

- ▶ **Electronic magnifiers**
- ▶ **Video magnifiers**
- ▶ **Camera systems**
- ▶ **Reading systems**
- ▶ **Accessories: XY-tables, height adjustable tables, etc.**

Product range





Reading systems

Reading systems help highly visually impaired or blind people to read text individually and independently. The systems convert text to speech and read them out with a clear voice. Volume and speed of the voice can be adjusted individually.

Reading systems capture printed text with a scanner and digitise it by using a text to speech e.g. OCR-function*. Large buttons in different shapes help the user to operate the system with ease. Most reading systems are designed for the mobile use. Some systems are equipped with a low vision function.

Find out more about our reading system in this brochure.

* OCR = optical character recognition

Facilitate everyday life with modern reading systems



VOXbox

The VOXbox is a new modular and multifunctional reading system based on an intuitive operational concept.

The system can be ordered with or without battery.

When needed, the reading system can be expanded with a low vision-module. By connecting the VOXbox scanner to a PC or Laptop, the unit develops into a combined reading and magnification device with easy navigation with the use of the PC mouse.

With the push of a button, documents can be captured easily, within a few seconds, by the powerful scanner. Offering the user high flexibility, the VOXbox will read aloud with a clear, precise, natural voice. Reading speed and sound volume can be adjusted to accommodate an individuals needs.

Designed for ease and mobility, the VOXbox carrying case is ergonomic and will transport the VOXbox safely and conveniently to any destination. Upon opening, the VOXbox is ready for immediate start up. The removable control unit guarantees convenient, smooth operation.

The modular reading system

"For further information please contact your Reinecker Vision partner."



VOX_{box} – your benefits at a glance

- high quality reading functionality with individual adjustment of reading speed and sound volume as well as navigation per paragraph, sentence, word and letter
- powerful scanner will adjust automatically to fit a standard letter size document, which can be easily positioned and integrated to the VOX_{box}
- intuitive operation concept: removable controller, and ease of start up, expandable to work with PC
- modern smart design with included ergonomically correct transport case for ease of use and quick start up



expandable



intuitive



mobile



VOX_{box} mains operation

OR

VOX_{box} including battery

includes 2 speakers,
headphone jack and
headphones

weight approx.
5 lbs (without battery)
6.3 lbs (with battery)
dimensions (W x H x D)
approx. 11.8" x 10.6" x 5.1"

navigation:
by paragraph, sentence, word,
letter

with optional low vision
module: different color
combinations, 6 different
display modes

included items and features:
control unit, cover for control
unit, shoulder strap, USB
interface (4 connections),
aid for easy positioning of the
documents, headphones, user
manual

10 languages
with many voices,
e.g. English, French, Spanish,
Italian, Polish, Russian,
German, ect.

magnification with
optional low vision module:
up to approx. 110x
(with 24" monitor)

available option:
low vision-module,
2in1-option (for the use of the
camera scanner with the VOX-
viva-software-CD
on a PC)



Reinecker USA LLC

145 River Rock Drive
Buffalo, NY 14207
USA

Tel. 716-348-3517

Fax 413-403-9777

Toll free: 866-733-2352

sales@reineckerusa.com

www.reineckerusa.com