

Videomatic

The classic video magnifier in a new look



We're here for you

We at Reinecker Vision are happy and proud to keep our traditional, independent, family owned standing since 1971 and always put our eye on intelligent and innovative product developments, all made in Germany.

With all of our product innovations that we export to more than 50 countries, we focus first on the demands of our low vision individuals.

Whether in the home, educational or professional environment we design products to help our customers reach their full potential.

Your benefits:

- more than 45 years of experience in low vision solutions
- short decision making progress
- individual support
- one face to the customer care

Reinecker Vision

Everything from one source

If your vision declines considerably at an older age or as a result of an eye disease, a wide range of assistive devices can help you manage your everyday life. We are happy to help you put together individual solutions to improve your daily life.

Our product range for your well-being

- ▶ **Electronic magnifiers: MANO**
- ▶ **Video magnifiers: mezzo, VEO, Videomatic**
- ▶ **Camera systems: topolino**
- ▶ **Reading systems: VOX**
- ▶ **Accessories: XY-tables, height adjustable tables, etc.**

Product range





Video magnifiers

Video magnifiers are reading aids mostly for stationary use. They capture text or pictures and reproduce it on a screen with a higher magnification. Video magnifiers offer a maximum of comfort and magnification.

Important features of video magnifiers are a high resolution and an optimized illumination with reduced reflections. Due to an auto focus the camera adjusts automatically to the optimum sharpness. To a certain extent objects and thick books will be displayed sharply despite of different heights.

In addition to full color display, video magnifiers offer a variety of semi color combinations in order to obtain a contrast enhancement. Text and objects can easily be moved underneath the camera with a reading table (XY-table).

You find more detailed information about the Videomatic on the following pages.

Reading comfortably with desktop CCTVs

Zubereitung

Das Fleisch waschen, trocken
2 Esslöffel Olivenöl in einem
ten. Möhren und Sellerie zuge
mian einrühren



Videomatic

In 1971, the first video magnifier manufactured by Reinecker Vision, was named Videomatic. Ever since it stands for comfort, innovation and quality.


We are now presenting Videomatic in a new design and combining proven quality features with contemporary design and innovative technology.

The high-resolution flat screen is mounted on an ergonomic swivel-tilt arm that helps the user to position it ergonomically according to individual needs. The operating panel is clearly structured with four large, tactile controls that make operation of Videomatic easy and intuitive.

Videomatic is equipped with a fullHD camera which magnifies texts of books or newspapers and other printed material such like letters. You can easily write or do cross words on the movable XY-table with friction brake. The three-level light adjustment further aids in adjusting Videomatic to your personal needs.

The contemporary video magnifier

The Focus camera – with technology

- fullHD camera with digital zoom and comfort autofocus that can be turned OFF
- magnification by default from approx. 4x to 42x (refers to factory settings, further magnification levels to be adjusted within the internal menu)
-  technology offers a high contrast, flicker free picture in HD+ quality. The innovative low vision camera has been edited optimally to reduce noise and blurring effects of the picture.



ergonomical



comfortable



intuitive

Videomatic – your benefits at a glance

- high resolution 20"-class LED flat screen (approx. 48 cm in diagonal)
- ergonomic and comfortable positioning of the height adjustable and inclinable monitor
- three level light adjustment
- individual positioning of text or pictures underneath the camera due to movable XY-table with friction brake
- ease of operation via four clearly structured and tactile controls
- overview mode: at the push of a button text will be zoomed out for better orientation at higher magnification levels
- line masking and shadowing offer a targeted orientation on screen
- solid, stable and light weight (approx. 12 kg/26.5 lbs)
- internal menu for further individual and customized adjustments



LED flat screen

20" class | approx. 49cm

adjustable in height,
inclinable



Focus camera: FullHD camera
with digital zoom

magnification approx. 4x to 42x
*

*refers to factory settings, further magnification
levels to be adjusted within the internal menu



true colors

semi colors:



individually selectable
semi color modes



scope of supply:
comfortable XY table,
line masking, shadowing,
overview mode,
dimming function,
individual menu



weight approx. 12 kg/26.5 lbs

dimensions approx. (W x H x D)
46.5 x 54 x 49 cm/
18.3" x 21.3" x 19.3"



options:
height adjustable table



Reinecker USA LLC

145 River Rock Drive
Buffalo, NY 14207
USA

Tel. 716-348-3517

Fax 413-403-9777

Toll free: 866-733-2352

sales@reineckerusa.com

www.reineckerusa.com