

# VEO UNO

Reimagined – in every detail



The VEO UNO video magnifier offers maximum reading and operating comfort in a modern, fresh design.

As a logical evolution of the proven VEO, the VEO UNO impresses with innovative technology and a well-thought-out design—developed and manufactured in Germany.

NELE Technology 2.0 minimizes reflections and shadows almost completely, even with high-gloss documents, thereby creating noticeably better image quality. The optimized XY table allows for more precise movement of the document, and the compact design frees up more space on your desk.

## Your benefits at a glance

- **Flexible magnification:** From approx. 1.9x (approx. 3.1x factory setting) up to over approx. 110x (approx. 44x factory setting).
- **Always in focus:** The optimized focus camera with an exceptionally wide depth of field eliminates the need to turn off autofocus while writing.
- **NELE Technology 2.0:** Virtually reflection-free image for noticeably better image quality.
- **Intuitive haptics:** High-contrast, easily accessible buttons for all settings.
- **Optimized XY table:** For even more precise and controlled work.
- **Compact build:** Modern, foldable design with a very small footprint.
- **Sustainable concept:** Highly durable, repair-friendly, equipped with a USB-C power adapter.
- **Cleverly packed:** Simplified system for quick and easy storage.

# VEO UNO

order code: VEOUNO241



LED flat screen  
24" class | approx. 60 cm (1920x1080 px)



Focus camera:  
FullHD camera with digital zoom  
comfortable autofocus with one-touch  
deactivation



from approx. 1.9x (factory setting:  
approx. 3.1x) up to over approx. 110x  
(factory setting: approx. 44x)



true colors,  
semi colors: **F F F F F**  
individually selectable semi color modes



weight approx. 10.5 kg  
dimensions (W x H x D) approx.  
54 x 54 x 54 cm,  
folded: approx. 22 cm high



Scope of supply:  
USB-C power adapter with cable,  
user manual



high-contrast controls  
Comfortable reading with an XY table



space-saving, foldable



ergonomically adjustable monitor

**Reinecker Vision GmbH**  
Head office Germany  
Sandwiesenstraße 19  
D-64665 Alsbach-Hähnlein

**Tel. +49 6257 9311-0**

Fax +49 6257 9311-8008  
international@reineckervision.de

**Your authorized distributor:**